


Specifications



FS-L1901D
19"



FS-L2401D
24"



FS-L2601D
26"



FS-L3201D
32"



FS-L4202D
42"

	FS-L1901D	FS-L2401D	FS-L2601D	FS-L3201D	FS-L4202D
TFT LCD Panel	19" LG Display	24" LG Display	26" LG Display	32" LG Display	42" LG Display
Panel Type	IPS				
Resolution	1280 x 1024	1920 x 1200		1920 x 1080	
Aspect Ratio	5:4	16:10		16:9	
Brightness (cd/m2)	450	400	500		
Contrast Ratio	800:1		1000:1	1300:1	
Response Time (Tr + Tf) ms	5	12	12	8	8
Viewing Angles (U/D/L/R)	89/89/89/89				
Image Scaler	Pixelworks				
Color Saturation	100% SMPTE 286M				
Dimensions (W x H x D) mm	442.3 x 338.2 x 76.2	580 x 386 x 87.5	618 x 412 x 99.5	795 x 485 x 100	1024.6 x 617.4 x 111.1
Enclosure Rating	IPX1				
Power Supply Type	External	Ext. Locking connector	Ext. Locking connector	Internal	Internal
Weight	14lbs. (6.4 kg)	20lbs. (9kg)	20lbs. (8.8kg)	31lbs (14kg)	66lbs. (30kg)
Operating Temp.	32 to 104° F (0 to 40° C)				
Available Inputs	(3 = Yes, N / A = Not Available)				
HD/SD-SDI	3	3	3	3	3
RGBS	3	3	3	3	3
RGBHV	3	3	3	3	3
	3	3	3	3	3
VGA	3	3	3	3	3
Component	3	3	3	3	3
S-Video	3	3	3	3	3
Composite	3	3	3	3	3
Fiber Optic DVI	3	3	3	3	3
RS-232C	3	3	3	3	3
Available Outputs	(3 = Yes, N / A = Not Available)				
HD/SD-SDI	3	3	3	3	3
DVI	3	3	3	3	3
Special Features	(3 = YES, N / A = Not Available, O = Optional)				
PIP/PBP	3	3	3	3	3
PAN/ZOOM	3	3	3	3	3
FREEZE	3	3	3	3	3
Touch Screen	3	3	3	3	3
GPIO	3	3	3	3	3
Custom Firmware	O	O	O	O	O
Remote		N/A		3	3
Certifications	(3 = YES, P = Pending)				
UL60601-1	P	3	3	3	3
EN60601-1-1	P	3	3	3	3
EN60601-1-2	P	3	3	3	3
CE, CCC, FCC-B	P	3	3	3	3
CAN/CSA	P	3	3	3	3
ICES-001	P	3	3	3	3
VCCI Class B	P	3	3	3	3
Medical Class 1	P	3	3	3	3

